



DNA GEL ELECTROPHORESIS SYSTEM

Rabrotech horizontal gel electrophoresis system provides the best quality in innovative design and manufacture, acquired during our existence. Each model incorporates the many design and safety features, either from within our in house product development team or laboratories worldwide. As a result, all our models are manufactured and finished to the highest specification, easy to use, and with safety always paramount.

Horizontal Gel units

Mini Horizontal	(Gel Size 6X8cm)
Mini Plus Horizontal	(Gel Size 10X12cm)
Midi Horizontal	(Gel Size 15X15cm)
Midi Plus Horizontal	(Gel Size 15X19cm)
Maxi Horizontal	(Gel Size 25X22cm)
Maxi Plus horizontal	(Gel Size 25X33cm)





Mini Horizontal DNA Gel Electrophoresis Unit RE001

Cat No RE001 Mini Horizontal DNA Gel Electrophoresis System

DNA Electrophoresis. Gel electrophoresis separates DNA fragments by size in a solid support medium (an Agarose gel.)

Minimum Sample	8 wells
Maximum sample	32wells
Gel Size	6x8cm
Buffer Volume	480ml
Comb slots distance	3.5cm
Voltage	75-95V

Products Code	Product	Size	Qty.
RE001-UT	Mini horizontal gel running tank 13.5 x 26 x 6cm	480ml	1
RE001-RT	Gel Casting Tray (Gel running tray)	6X8 cm	1
RE001-CT	Gel caster tray with silicon mat	8.5x14cm	1
RE001-PLUS-W	0.2mm thick platinum wire	12cm	1pair
RE15C	Power leads with 4mm output connectors	1.5 meter	1pair
RE001-C8-1	8Well comb, teeth size 4.5mm,	1mm	2
RE001- C16-1	16 Well comb, teeth size 2.2mm	1mm	2
RE001-C8-1.5	8well comb, teeth size 4.5mm	1.5mm	2
RE001-C16-1.5	16 Well comb ,2.2mm	1.5mm	2

Optional Accessories:

Products Code	Product	Size	Qty.
RE001-RT	Gel Casting Tray (Gel running tray)	6X8 cm	1
RE001-CT	Gel caster for 6x8cm tray	8.5X13cm	1
RE001-PLUS-WR	0.2mm thick platinum wire Red	12cm	1
RE001-PLUS-WB	0.2mm thick platinum wire Black	12cm	1



RABROTECHTM
Science for a better life

Products Code	Product	Size	Qty.
RE15C	Power leads with 4mm output connectors	1.5 meter	1pair
RE001-C8-1	8Well comb, teeth size 4.5mm,	1mm	1
RE001-C12-1	12 Well comb, teeth size 2.5mm	1mm	1
RE001- C16-1	16 Well comb, teeth size 2.2mm	1mm	1
RE001-C8-1.5	8Well comb, teeth size 4.5mm	1.5mm	1
RE001-C12-1.5	12 Well comb, teeth size 2.5mm	1.5mm	1
RE001-C16-1.5	16 Well comb, teeth size 2.2mm	1.5mm	1
RE001-C8-2	8Well comb, teeth size 4.5mm	2mm	1
RE001-C12-2	12 Well comb, teeth size 2.5mm	2mm	1
RE001-C16-2	16 Well comb, teeth size 2.2mm	2mm	1



RABROTECHTM
Science for a better life

B-352 Rohini Sector-1 Avantika Delhi-110085

Mail Id: info@rabrotech.com

Ph No: 9999914074

WWW.rabrotech.com