



DNA GEL ELECTROPHORESIS SYSTEM

Rabrotech horizontal gel electrophoresis system provides the best quality in innovative design and manufacture, acquired during our existence. Each model incorporates the many design and safety features, either from within our in house product development team or laboratories worldwide. As a result, all our models are manufactured and finished to the highest specification, easy to use, and with safety always paramount.

Horizontal Gel units

Mini Horizontal	(Gel Size 6X8cm)
Mini Plus Horizontal	(Gel Size 10X12cm)
Midi Horizontal	(Gel Size 15X15cm)
Midi Plus Horizontal	(Gel Size 15X19cm)
Maxi Horizontal	(Gel Size 25X22cm)
Maxi Plus horizontal	(Gel Size 25X33cm)





Mini Plus Horizontal DNA Gel Electrophoresis Unit RE001PLUS

Cat No RE001PLUS Mini Plus Horizontal DNA Gel Electrophoresis System

DNA Electrophoresis. Gel electrophoresis separates DNA fragments by size in a solid support medium (an Agarose gel.)

Minimum Sample	10 wells
Maximum sample	40wells
Gel Size	10x12cm
Buffer Volume	480ml
Comb slots distance	3.5cm
Voltage	75-120V

Products Code	Product	Size	Qty.
RE001- PLUS-UT	Mini Plus horizontal gel running tank 13.5 x 26 x 6cm	480ml	1
RE001-PLUS-RT	Gel Casting Tray (Gel running tray) with two slots	10X12cm	1
RE001- PLUS -CT	Gel caster tray with silicon mat	12x17cm	1
RE001- PLUS -W	0.2mm thick platinum wire	12cm	1pair
RE15C	Power leads with 4mm output connectors	1.5 meter	1pair
RE001-PLUS-C10-1	10Well comb, teeth size 8mm/2.5mm teeth space	1mm	2
RE001- PLUS- C16-1	16 Well comb, teeth size 3.5mm/2mm teeth space	1mm	2
RE001- PLUS-C10-1.5	10well comb, teeth size 8mm/2.5 teeth space	1.5mm	2
RE001- PLUS C16-1.5	16 Well comb ,3.5mm//2mm teeth space	1.5mm	2

Optional Accessories:

Products Code	Product	Size	Qty.
RE001-PLUS-RT	Gel Casting Tray (Gel running tray) with two slots	10x12cm	1
RE001-PLUS-CT	Gel caster tray with silicon mat	12x17cm	1

Products Code	Product	Size	Qty.
RE001- PLUS -WB	0.2mm thick platinum wire Black	12cm	1
RE001- PLUS -WR	0.2mm thick platinum wire Red	12cm	1
RE15C	Power leads with 4mm output connectors	1.5 meter	1pair
RE001-PLUS-C10-1	10Well comb, teeth size 8mm/2.5mm teeth space	1mm	1
RE001- PLUS- C16-1	16Well comb, teeth size 3.5mm/2mm teeth space	1mm	1
RE001-PLUS-C20-1	20Well comb, teeth size 3mm/1.5mm teeth space	1mm	1
RE001-PLUS-C10-1.5	10Well comb, teeth size 8mm/2.5mm teeth space	1.5mm	1
RE001-PLUS- C16-1.5	16Well comb, teeth size 3.5mm/2mm teeth space	1.5mm	1
RE001-PLUS-C20-1.5	20Well comb, teeth size 3mm/1.5mm teeth space	1.5mm	1