




DNA GEL ELECTROPHORESIS SYSTEM

Rabrotech horizontal gel electrophoresis system provides the best quality in innovative design and manufacture, acquired during our existence. Each model incorporates the many design and safety features, either from within our in house product development team or laboratories worldwide. As a result, all our models are manufactured and finished to the highest specification, easy to use, and with safety always paramount.

Horizontal Gel units

Mini Horizontal	(Gel Size 6X8cm)
Mini Plus Horizontal	(Gel Size 10X12cm)
Midi Horizontal	(Gel Size 15X15cm)
Midi Plus Horizontal	(Gel Size 15X19cm)
Maxi Horizontal	(Gel Size 25X22cm)
Maxi Plus horizontal	(Gel Size 25X33cm)





Maxi Plus Horizontal DNA Gel Electrophoresis Unit RE003PLUS

Cat No RE003-PLUS MAXI PLUS Horizontal DNA Gel Electrophoresis System

DNA Electrophoresis. Gel electrophoresis separates DNA fragments by size in a solid support medium (an Agarose gel.)

Minimum Sample	26 wells
Maximum sample	312wells
Gel Size	25x33cm
Buffer Volume	3Ltr.
Comb slots distance	2cm
Voltage	150-200V

Products Code	Product	Size	Qty.
RE003-PLUS-UT	Maxi Plus horizontal gel running tank 29x53x8cm	3Ltr.	1
RE003- PLUS-RT	Gel Casting Tray (Gel running tray) with 12 slots	25x33cm	1
RE003- PLUS-CT	Gel caster tray with silicon mat	28x40cm	1
RE003- PLUS-W	0.3mm thick platinum wire	27cm	1pair
RE15C	Power leads with 4mm output connectors	1.5 meter	1pair
RE003-PLUS-C26-1	26 Well comb, teeth size 7mm/2mm teeth space	1mm	4
RE003- PLUS-C26-1.5	26 Well comb, teeth size 7mm/2mm teeth space	1.5mm	4

Optional Accessories:

Products Code	Product	Size	Qty.
RE003- PLUS-RT	Gel Casting Tray (Gel running tray) with 12 slots	25x33cm	1
RE003- PLUS-CT	Gel caster tray with silicon mat	28x40cm	1
RE003- PLUS-WB	0.3mm thick platinum wire Black	27cm	1
RE003- PLUS-WR	0.3mm thick platinum wire Red	27cm	1
RE15C	Power leads with 4mm output connectors	1.5 meter	1pair



RABROTECHTM
Science for a better life

Products Code	Product	Size	Qty.
RE003-PLUS-C26-1	26 Well comb, teeth size 7mm/2mm teeth space	1mm	1
RE003-PLUS-C26-1.5	26 Well comb, teeth size 7mm/2mm teeth space	1.5mm	1
RE003-PLUS-C52-1	52 Well comb, teeth size 3mm/1.5mm teeth space	1mm	1
RE003-PLUS-C52-1.5	52 Well comb, teeth size 3mm/1.5mm teeth space	1.5mm	1



RABROTECHTM
Science for a better life

B-352 Rohini Sector-1 Avantika Delhi-110085
Mail Id: info@rabrotech.com
Ph No: 9999914074
WWW.rabrotech.com